

PRESSALIT.



Shower and changing tables

Pressalit shower and changing tables are designed and manufactured to provide a reliable foundation for the care of children and adults. Our shower and changing tables are fashioned to aid the comfort of the caregiver and care receiver.

Questions or in need of advice? Call us today at

317-494-2889

Reliable solution Designed with people in focus

According to the World Report on Disability, millions and millions of people around the world are living with some form of physical disability. This makes them the world's largest minority group. Yet their needs are frequently overlooked in schools, hospitals, and other public spaces, mainly due to inaccessible or inadequate restroom facilities.

Dignity

For people with physical or mental difficulties, a simple trip may be impossible – or humiliating and undignified – because the restroom facilities fail to accommodate their needs.

Due to the lack of facilities and assistive equipment in restrooms, many intentionally stop consuming fluids before leaving home to avoid inadequate public restrooms. Others are forced to lie down on unhygienic restroom floors or in the back of their cars for care or changing. This is why accessible restrooms with shower and changing tables are essential in all public places.

Reliable safety for users and caregivers

Optimal working conditions are vital for providing care to people living with disabilities. Caregivers and family members are often subjected to dangerous working conditions when providing care in a facility not properly equipped to meet their needs. This can cause back, shoulder, and other injuries, putting the safety of the user and caregiver at risk.

Pressalit height adjustable shower and changing tables ensure optimal working conditions for caregivers and family members. Safety rails provide privacy and comfort for the user. The height adjustable feature makes lifting to and from the table easier and helps to avoid uncomfortable and dangerous working positions.



Easy to clean

Maintaining excellent hygiene in restrooms require fixtures that are quick and easy to clean. Pressalit shower and changing tables are designed and built with a strong focus on healthy living. This makes the shower and changing tables ideal for use in public institutions such as hospitals and schools.

Flexibility and quality

The right restroom facilities contribute to an increased sense of independence for people with physical disabilities. The higher degree of self-reliance also improves the user's overall quality of life.

Pressalit offers a lot of features to adapt shower and changing tables to match a given space and meet individual user needs. All Pressalit shower and changing tables come with a 3-year guarantee – our assurance of a stable, durable, and high quality solution.

EXPLORE OUR PRODUCT RANGE ONLINE

READ MORE ABOUT HOW WE TEST THE QUALITY AND SAFETY IN OUR PRODUCTS ONLINE



Technology like this not only improves quality of life for users – it also contributes to a healthier working environment for family and caregivers"

Holger Søe Plougsgaard, Head of Product Development at Pressalit

PITTSBURGH INTERNATIONAL AIRPORT Becomes all accessible

As the first airport in the world, Pittsburg International Airport has introduced a dedicated facility to help reduce stress for travelers with sensory sensitivity. Here, those with autism or similar issues and their families can de-escalate and relax prior to take-off and after landing. The Pittsburgh International Airport chose Pressalit products to improve accessibility in the terminal's family restroom. They have installed a Pressalit height adjustable sink and height adjustable adult-sized changing table. This makes it possible for families to travel comfortably, regardless of physical and mental ability.



"

It is so important to make sure that I am not injuring my family members or my husband when they are trying to make sure that my needs are taken care of. The safety of both me and the care provider is vital and that is one reason the height adjustable tables are so important."

Dr. Josie Badger, diagnosed with congenital myasthenic syndrome



Facilities to meet the needs of the community

Nothing was left to chance when the Pittsburgh International Airport began the process of making the airport more accessible. This was a strong commitment to establish facilities that meet the needs of the entire community. The airport consulted with advocacy groups, caregivers of children with neurodevelopmental challenges, and other experts to form a winning plan.

This consultative process collected experiences and suggestions for the design of the sensory room and the family restroom. It became clear that one of the most important items for users is the availability of an adult-sized changing table – and that it is height adjustable. This improves the safety of use and eases the transfer from wheelchair to changing table. Travelers with disabilities would now have the dignity and privacy they need when visiting the airport.


A better travel experience

Many families impacted by a disability do not fly long distance – in fact many do not fly at all. One reason for this is the lack of accessible facilities in family restrooms. The installation of a Pressalit adult-sized changing table and height adjustable sink is a major step to improve restroom accessibility. The impact is substantial. More people than expected use the changing table and the height adjustable sink, which improves the airport experience for visitors. Travelers have less stress and are more comfortable, which makes the flight more relaxing for everyone.

GET INSPIRATION ONLINE

SEE MORE REFERENCES ONLINE





Mobile shower and change trolley A flexible, safe, and easy-to-clean solution

Many settings need mobile shower and changing solutions that can be moved from place to place and at the same time accommodate the range of needs arising for users, caregivers, and family members. Pressalit's change trollleys can easily be adapted to a specific care situation and create space for flexibility in a safe and easy-to-clean care environment.

Because the mobile shower change trolley and mobile change trolley were both designed for the purpose of creating the best possible solution for users and caregivers, every detail and function on the trolleys have been carefully planned and tested by the designer in collaboration with professional care personnel in the development phase.

A comfortable care situation requires secure solutions

Both the shower trolley and change trolley are equipped with safety rails along the sides, while the shower trolley also features a headboard and footboard. For individuals with sensory limitations, the firm, continuous safety rails help define the space visually. This not only helps prevent fall accidents; their robust extension and secure construction also serve to make the entire care situation safe and comfortable for the user. Furthermore, this helps reducing the risk of muscle tone disorders.



Easy operation creates security for user and caregiver

For caregivers and family members providing care around the trolley, it is important that the trolley is both user-friendly and safe. This is the key element to creating an essential feeling of security. The rails, providing safety and increasing the user's sense of psychological well-being, are easily operated with one hand. This makes the caregiver's job easier. He or she can maintain uninterrupted contact with the user for increased security from any position around the trolley. Because the mobile shower and change trolley can be moved around and the safety rails can be lowered, the caregiver can easily ensure optimal working conditions on each side.

The user can lie on the trolley comfortably and securely while being moved, thereby reducing the number of shifts. The large diameter wheels ensure easy, comfortable transportation on a range of surfaces, and the central lock system and non-slip foot pedals make it easy to lock all wheels at the same time.

The mattresses, which come separately, add to the feeling of comfort for the user when lying on the trolley, and at the same time they contribute to making the care situation easier for the caregiver. The lying surface being embraced by the folded mattress on the shower trolley creates a basin, which is easy to drain and quick to dry. This allows the user to remain lying on the trolley throughout the entire drying and dressing process, without water forming puddles on the mattress. Additionally, the shower as well as the change trolley mattress is easy to clean, thereby effectively maintaining optimal hygiene.

SEE PRODUCT SUMMARY ONLINE

READ MORE ABOUT THE MOBILE SHOWER/CHANGE TROLLEY ONLINE




• Working together with trained care personnel, my own body was lifted, bathed, dried and dressed. This formed the basis of my approach to design optimal and comfortable solutions for both caregivers and users.

Anne Mette Tygesen, Product Designer, Pressalit A/S

Overview Shower and changing tables





Product no. R8592318301

SCT 3000 Shower and change table with safety rail, electrically height adjustable Outer measurements: 41.6" x 75.2"

- Also available in 55.5"
- Maximum load: 440 lbs
- Also available with sapphire blue pads

Product no. R8435313299

SCT 3000 Shower and change table, with safety rail, fixed height. Outer measurements: 40.6" x 55.5"

- Also available in 75.2"
- Maximum load: 440 lbs
- Also available with graphite grey pads

Product no. R8592418000

SCT 2100 Shower and change table with safety rail, canvas cover, electrically height adjustable Outer measurements: 35.4" x 70.9"

- Also available with fine-mesh net
- Maximum load: 385 lbs

Product no. R8402000-93

SCT 1000 Shower and change table with safety rail, electrically height adjustable Outer measurements: 30.1" x 55.1"

- Also available in 70.9"
- Maximum load: 330 lbs



PRESSALIT.







Product no. R8595572000

CT 4100 Change table, electrically foldable, electrically height adjustable, with buttons placed in table front. Including safety rail and safety belt. Outer measurements: 40.2" x 72"

- Also available with manual folding mechanism
- Mattress sold separately

Product no. R9464721000

MSCT 1 Shower change trolley, electrically tilt and height adjustable, central lock system, storage drawer, four one-hand operated fold-down sides, battery-powered. Outer measurements: 32.3" x 86.2"

- Also available in 32.3" x 70.5"
- Shower mattress sold separately

Product no. R9462721000

MCT 2 Change trolley, electrically height adjustable, central lock system, storage drawer, two one-hand fold-down siderails, battery-powered. Outer measurements: 32.3" x 86.2"

- Also available in 32.3" x 70.5"
- Mattress sold separately

A shower and change table from Pressalit

- Provides flexibility to meet challenges through different stages of illness and for a wide variety of users.
- Is easy to clean, enabling good hygiene.
- Ensures optimal working conditions for family members, caregivers, and other personnel.
- Can be adapted to a given room or specific user.
- Provides peace of mind, comfort and safety for users with physical and cognitive disabilities.
- Makes the care situation dignified and comfortable for the user.
- Is an ideal solution for private homes, public restrooms, institutions, schools, hospitals, airports and hotels.



Pressalit A/S Pressalitvej 1 8680 Ry Denmark T: +45 8788 8788 F: +45 8788 8789 pressalit@pressalit.com

Pressalit Inc.

30 Westgate Parkway Ste 116 Asheville, NC 28806 USA T: +1 202 271 0166 us@pressalit.com

Pressalit.com

02.2020

PRESSALIT SUPPORTS UNITED NATIONS GLOBAL GOALS



PREMIUM BATHROOM SOLUTIONS SINCE 1954

Pressalit A/S is a Danish privately owned company with production facilities in Ry, Denmark. Since 1954, we have been designing and producing toilet seats of the highest quality. And since 1975, Pressalit has also been developing accessible bathrooms and height adjustable kitchen solutions for people with reduced mobility. This has made us one of the world's leading manufacturers of premium bathroom solutions. High quality, flawless design and reliable functionality are the core features of our wide range of products.



Reservations are made for typographic errors, technical changes and model deviations.